

DURATWIX® UNI LINGO

- ▶ Tracheostomy tube with wide inner lumen to provide maximum volume of air at optimized flow
- ▶ Cannula with 15 mm connector (UNI)
- ▶ Soft double-swivel ROTATWIX® neckflange, rotating in 360° for enhanced neck mobility
- ▶ With ASPILO® sieve-fenestration - roof-tile shaped (horizontal)
- ▶ Medical-grade polyurethane, soft, atraumatic, thermosensible



SHORT


Order information:

DURATWIX® UNI LINGO SIZE 8-10



REF 19271

DURATWIX® UNI LINGO short SIZE 7-9



REF 19272

When ordering, please state the size after the item number!

Size	O.D. neck flange mm	O.D. Tip mm	I.D. mm	Length mm standard	REF standard	Length mm short	REF short
7	11,9	10,8	7,0	-	-	72,0	19272-07
8	12,5	11,4	8,0	85,0	19271-08	75,0	19272-08
9	13,7	12,5	9,0	88,0	19271-09	78,0	19272-09
10	15,0	13,8	10,0	91,0	19271-10	-	-

O.D. neck flange = Outer diameter behind the neck flange · O.D. Tip = Outer diameter of the outer tube at the tip of the cannula · I.D. = Inner diameter at the tip of the cannula · Length = Length over the outer curve