

DURATWIX® UNI CUFF

- ▶ Tracheostomy tube with wide inner lumen to provide maximum volume of air at optimized flow
- ▶ Tracheostomy tube with 15mm connector (UNI), without inner cannula
- ▶ Soft double-swivel ROTATWIX® neckflange, rotating in 360° for enhanced neck mobility
- ▶ Medical-grade polyurethane, soft, atraumatic, thermosensible
- ▶ MUCOPROTECT® cuff: very thin walls for proper seal, reduced surface contact with tracheal Mucosa to minimize pressure area
- ▶ Set includes broad and adjustable tube holder, obturator and product ID card
- ▶ Size 7–9 available in two lengths



Order information:

DURATWIX® UNI CUFF, SIZE 7-9



REF 19601

DURATWIX® UNI CUFF short, SIZE 7-9



REF 19602

When ordering, please state the size after the item number!

Size	O.D. neck flange mm	O.D. Tip mm	I.D. mm	Length mm standard	REF standard	Length mm short	REF short
7	11.9	10.8	7.0	82.0	19601-07	72.0	19602-07
8	12.5	11.4	8.0	85.0	19601-08	75.0	19602-08
9	13.7	12.5	9.0	88.0	19601-09	78.0	19602-09

O.D. Tip = Outer diameter of the outer tube at the tip of the cannula · O.D. neck flange = Outer diameter behind the neck flange · I.D. = Inner diameter at the tip of the cannula · Length = Length over the outer curve